

# Alpha Support Belt


Light  
Easy to use  
High SWL



Support Belt



Support Belt Soft



The Support Belts are a part of a product range with a wide range of possible applications. This product is used in particular when raising and/or transferring clients in home care and other health care facilities and hospitals. It is a safe transfer aid, which allows for a pleasant and comfortable transfer for the client and being back-friendly and straightforward for the caregiver.


The interior of the Support Belt is lined with a non-slip material, so it clings to clothing, and vertical and horizontal grips make the handling easier. The robust outer material is made of nylon and prolongs the service life of the belt. The model that is called SOFT has a soft, velvet-like, inner material.

## Technical specification

Product Name	Art. nr.	Dimensions	Material	Max. user weight
Support Belt Extra Small	PM-6012	750 x 120 x 7 mm / 29.5" x 4.7" x 0,2"	Polyamide, latexfree	200 kg / 400 lb
Support Belt Soft Extra Small	PM-6012 Soft			
Support Belt Small	PM-6013	750 x 120 x 7 mm / 29.5" x 4.7" x 0,2"	Polyamide, latexfree	200 kg / 400 lb
Support Belt Soft Small	PM-6013 Soft			
Support Belt Medium	PM-6014	1170 x 120 x 7 mm / 43.9" x 4.7" x 0,2"	Polyamide, latexfree	200 kg / 400 lb
Support Belt Soft Medium	PM-6014 Soft			
Support Belt Large	PM-6015	1170 x 120 x 7 mm / 43.9" x 4.7" x 0,2"	Polyamide, latexfree	200 kg / 400 lb
Support Belt Soft Large	PM-6015 Soft			
Support Belt Extra Large	PM-6016	1850 x 120 x 7 mm / 72.8" x 4.7" x 0,2"	Polyamide, latexfree	200 kg / 400 lb
Support Belt Soft Extra Large	PM-6016 Soft			
Looped Belt (accessory)	PM-3010			

**PETERMANN**  
part of the Human Care Group

These products are distributed by  
Human Care AB for Petermann GmbH.



= 200 kg / 440lb



Visit [www.humancaregroup.com](http://www.humancaregroup.com) for the latest product information and contact details to national offices and partners.