


Sliding Board

Smooth
Flexible
Sideward
transfer



PM-5024 / PM-5025



BM-5000



The Blue Motion Sliding Board is flexible and has a very smooth surface in order to minimize resistance due to friction. It is used to transfer clients, e. g. from the edge of a bed to a wheelchair, toilet or car seat.

The client does not have to be able to stand up, but can make the change in position while sitting, either with or without assistance.



Technical specification

Product Name	Art nr	Dimensions	Material
Blue Motion Sliding Board	BM-5000	74 x 24 cm / 29,1 x 9,4"	Polyethylene
Alpha Sliding Board	PM-5024	45 x 20 / 17,7 x 7,8"	Polyethylene
Alpha Sliding Board	PM-5025	60 x 25 / 23,6 x 9,8"	Polyethylene



 = 200 kg / 440lb

PETERMANN
part of the Human Care Group

These products are distributed by
Human Care AB for Petermann GmbH.



Visit www.humancaregroup.com for the latest product information and contact details to national offices and partners.